
CLEAR VISION

SOUND STRATEGIES

SOLID PERFORMANCE

# Versa Studio 操作指南

Manual of Versa Studio

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**Princeton  
Applied  
Research**

<http://www.ametek.com>

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**The VersaStudio makes "Echem Easy!"**

More Options. More Solutions.

<http://www.par-solartron.com.cn>

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1

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


*Two-Electrode Connection*

- Batteries
- Capacitors
- Fuel Cells
- Sensors
- Resistors

*Three-Electrode Connection*

- Aqueous Electrochemistry
- Corrosion
- EIS


*Four-Electrode Connection*

- H-Cell Setup
- Liquid –liquid interface

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
## 4 Versa Studio软件操作

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The screenshot displays the Versa Studio software interface with several key components labeled:

- 菜单栏 (Menu Bar):** Located at the top, containing options like Experiment, Data, View, Tools, Security, Window, and Help.
- 工具栏 (Toolbar):** Located below the menu bar, containing various icons for file operations and data management.
- 数据图形显示窗口 (Data Graph Display Window):** Contains two plots:
  - Zim (ohms) vs Zre (ohms):** A Nyquist plot showing a semi-circular arc in the upper half-plane.
  - |Z| (ohms) vs Frequency (Hz):** A Bode magnitude plot showing a peak in impedance at a specific frequency.
- 数据列表窗口 (Data List Window):** A table titled "Data View - 31 Points (All)" showing experimental data points.
 

Potential (V)	Current (A)	Elapsed Ti...
3.586 mV	-13.62 $\mu$ A	0.811
3.559 mV	-13.293 $\mu$ A	1.231
3.633 mV	-19.66 $\mu$ A	1.922
3.628 mV	-19.824 $\mu$ A	2.353
3.606 mV	-19.248 $\mu$ A	2.714
3.625 mV	-19.608 $\mu$ A	3.0849999
3.602 mV	-22.292 $\mu$ A	3.485
3.589 mV	-21.957 $\mu$ A	3.8649999
3.605 mV	-22.019 $\mu$ A	4.2149999
3.561 mV	-25.655 $\mu$ A	4.6159999
3.462 mV	-17.449 $\mu$ A	5.0169999
3.499 mV	-17.883 $\mu$ A	5.4179999
3.473 mV	-18.242 $\mu$ A	5.8189999
3.474 mV	-18.985 $\mu$ A	6.2199999
3.455 mV	-16.811 $\mu$ A	6.6209999
3.453 mV	-17.202 $\mu$ A	7.0219999
3.494 mV	-17.65 $\mu$ A	7.4229999
3.49 mV	-17.273 $\mu$ A	7.8239999
3.456 mV	-16.622 $\mu$ A	8.2249999
3.504 mV	-18.433 $\mu$ A	8.6259999
3.474 mV	-17.305 $\mu$ A	9.0269999
3.472 mV	-17.613 $\mu$ A	9.4279999
3.49 mV	-18.146 $\mu$ A	9.8289999
3.47 mV	-17.375 $\mu$ A	10.2299999
3.465 mV	-17.846 $\mu$ A	10.6309999
3.47 mV	-17.523 $\mu$ A	11.0319999
3.484 mV	-17.627 $\mu$ A	11.4329999
3.472 mV	-17.411 $\mu$ A	11.8339999
3.466 mV	-18.003 $\mu$ A	12.2349999
3.456 mV	-17.072 $\mu$ A	12.6359999
- 实验参数设置窗口 (Experiment Properties Window):** Located at the bottom, containing sections for:
  - Actions to be Performed:** Potentiostatic EIS.
  - Properties for Potentiostatic EIS:**
    - AC Properties: Start Frequency (10000 Hz), End Frequency (1 Hz), Amplitude (10 mV RMS).
    - DC Properties: Step or Scan (Step), Potential (0 V), Versus (vs OC).
    - Scan Properties: Point Spacing (Logarithmic), Number of Points (30), Points Per Decade (10), Data Quality (1), Measurement Delay (0 s).
  - Frequency List (Hz):** A list of frequencies ranging from 10000.000000 to 894.407700.
- 正常连接时显示工作站型号及其序列号 (Instrument: None):** Located at the bottom left of the interface.

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The screenshot shows the Versa Studio software interface. The 'Reference Electrodes' dialog box is open, displaying a list of electrodes and their voltages. The 'Add / Edit Reference Electrode' dialog box is also open, with red arrows pointing to the 'Electrode Name' and 'Voltage' input fields. The 'Tools' menu is open, highlighting 'Reference Electrode List...'. The 'Experiment Properties' window shows 'Linear Scan Voltammetry' as the action to be performed.

Electrode Name	Voltage (V)
Ag, AgCl / NaCl (sat'd)	0.194
Ag, AgCl / KCl (3.5M)	0.205
Ag / Ag <sup>+</sup> (0.1M) in acetonitrile	0.8
NHE Normal Hydrogen Electrode	0
SCE Saturated Calomel (sat'd KCl) (unspecified)	0.242
Ag, AgCl / KCl (sat'd)	0.197
SCE Saturated Calomel (sat'd NaCl)	0.236

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The screenshot shows the Versa Studio software interface with the Options dialog box open. The dialog box has two tabs: General and Instrument. The Instrument tab is selected, and several sections are highlighted with orange boxes and Chinese annotations:

- Current Polarity Convention:** Two radio buttons are shown. The second one, "Define cathodic current as negative", is selected. An annotation reads "电流极性定义, 国内标准" (Current polarity definition, domestic standard).
- Corrosion Properties:** Three input fields are shown: Cathodic Beta Constant (100 mV), Anodic Beta Constant (100 mV), and Corrosion Rate Units (mV). An annotation reads "腐蚀参数, 建议使用默认值" (Corrosion parameters, recommended to use default values).
- General Settings:** A group of checkboxes is shown, including "Enable Periodic Data Storage", "Restore Last Experiment", "Cell to Internal at Experiment End", "Lock Experiment Properties to Data", "Automatically Save Data File", "USB Compatibility Mode", and "Display Advanced Experiment Properties". An annotation reads "一般设置, 建议使用默认选项" (General settings, recommended to use default options).
- Language:** A dropdown menu is shown with "Chinese (Simplified, PRC)" selected. An annotation reads "界面语言可选择简体中文" (Interface language can be selected as Simplified Chinese).

The background shows the main software window with a graph of Current (A) vs Potential (V) and a menu where "Options..." is highlighted in red. The "Experiment Properties" panel at the bottom left shows "Linear Scan Voltammetry" as the action to be performed.

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Add Instrument to Group  
Select Instrument  
Locate Instrument  
Rename Instrument...  
Assign Setup File...  
Remove Setup File

1111111	1111113	1111118	1111110
E: 521.40 $\mu$ V	E: 61.34 $\mu$ V	E: 306.70 $\mu$ V	E: 214.69 $\mu$ V
I: 1.11 nA	I: 343.51 pA	I: 79.74 pA	I: 36.80 pA
Voltammetry Checkout	No Setup File	Voltammetry Checkout	No Setup File

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

Click **Common**

Select reference electrode and set working electrode parameters

Can choose whether to re-measure OCP before each experiment, especially for batch measurements of batteries

Safety limit, exceeding all limits will stop the experiment

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Properties	Value
Reference Electrode	(unspecified) (0 Volts)
Working Electrode Type	Unspecified
RDE Speed (Volts)	0
Working Electrode Area (cm <sup>2</sup> )	1
Density (g/ml)	0
Equivalent Weight (g)	0
Mass (g)	0
Re-measure OC Per Action	No
Measured Open Circuit	0

Limits	Direction	Value
None	≤	0
None	≤	0

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**Click the corresponding test method**

**Click "Advanced", for external testing Cell to Use must select External**

**Click "Help", you can view the specific meaning of the parameters of the test method**

**Set the specific parameters of the measurement method**

**For any 4-electrode measurement system, Differential must be selected. Current above 200mA, Differential can eliminate the error brought by the working electrode wiring resistance, with high accuracy. While Single Ended has higher stability, suitable for fluctuating systems, such as corrosion, high-capacity electrochemical cells and EIS measurement, etc.**

Step Properties	Value	Versus
Initial Potential (V)	0	vs Re
Final Potential (V)	1	vs Re

Scan Properties	Value
Scan Rate (V/s)	0.005
Total Points	1000

Limits	Direction	Value
None	≤	0
None	≤	0

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External
	External
	Internal

Graph: Voltage vs Time. Initial Potential (V) is marked at the start of the scan. Final Potential (V) is marked at the end. The scan rate is given as  $\text{Scan Rate (V/s)} = \frac{dE}{dt}$ . The Open Circuit (OC) potential is also indicated.

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Experiment Properties

Actions to be Performed:

- Common
  - Loop #1
    - Constant Current
    - Potentiostatic EIS
  - Loop #2
    - Cyclic Voltammetry
    - Chronoamperometry
    - Linear Scan Voltammetry

Properties for Linear Scan

Step Properties	Value	Versus
Initial Potential (V)	0	vs Ref
Final Potential (V)	1	vs Ref

Scan Properties	Value
Scan Rate (V/s)	0.005
Total Points	1000

Select an Action

Actions | Advanced Actions

Pre Experiment Actions	Sequence Actions
Condition	Loop
Deposition	Time Delay
Equilibration	Message Prompt
Purge	Measure Open Circuit
iR Determination	Auxiliary Interface
	Run External Application
	DAC Output Control
	EMail
	Auto Current Range Setup
	Display Message

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The screenshot displays the Versa Studio software interface. The main window shows two graphs: a Nyquist plot (Z<sub>im</sub> vs Z<sub>re</sub>) and a Bode magnitude plot (|Z| vs Frequency). The 'Add Graph' dialog box is open, showing a grid of available graph templates. The 'Nyquist' template is highlighted.

**Available Graph Templates:**

- Admittance
- Aux Chn 0 vs t
- Aux Chn 1 vs t
- Aux Chn 2 vs t
- Aux Chn 3 vs t
- Bode Phase
- Bode Z<sub>im</sub>
- Bode Z<sub>mag</sub>
- Bode Z<sub>re</sub>
- Capacitance
- Capacity vs t
- Corrosion Rate vs t
- Delta I (F-R) vs E
- Delta I (R-F) vs E
- E vs Energy
- E vs I
- E vs Log(I)
- E vs t
- Foward I vs E
- Frequency vs t
- I vs E
- I vs t
- Mott-Sch...
- Nyquist
- Power vs Energy
- Q vs t
- Resistance vs t
- Reverse I vs E
- Sync ADC Input vs t

**Experiment Properties:**

Common		AC Properties		DC Properties	
Start Frequency (Hz)	10000	Value		Value	Versus
End Frequency (Hz)	1			Step or Scan	Step
Amplitude (mV RMS)	10			Potential (V)	0 vs OC
				Final Potential (V)	0 vs OC
Scan Properties		Instrument Properties		Cell Properties	
Point Spacing	Logarithmic	Value		Value	
Number of Points	30			Current Range	Auto
Points Per Decade	10			Electrometer Mode	Single Ended
Data Quality	1			Bandwidth Limit	Auto
Measurement Delay (s)	0			LCI Bandwidth Limit	Auto
Scan Rate (V/s)	005			Cell Properties	Value
				Leave Cell ON	No
				Cell to Use	External

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










The screenshot shows the 'View' menu with 'Show Data View' selected. The 'Data View' window displays a table of experimental data for 50,000 points across all segments.


Potential (V)	Current (A)	Elapsed Tim...
999.249 mV	999.126 μA	1
999.249 mV	999.126 μA	2
999.249 mV	999.126 μA	3
999.249 mV	999.126 μA	4
999.249 mV	999.126 μA	5
999.249 mV	999.126 μA	6
999.249 mV	999.126 μA	7
999.249 mV	999.126 μA	8
999.249 mV	999.126 μA	9
999.249 mV	999.126 μA	10
999.249 mV	999.126 μA	11
999.249 mV	999.126 μA	12
999.249 mV	999.126 μA	13
999.249 mV	999.065 μA	14
999.249 mV	999.126 μA	15
999.249 mV	999.126 μA	16

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**Data Copy**

Columns to Copy to Clipboard:

Enable	Description
<input checked="" type="checkbox"/>	Potential (V)
<input checked="" type="checkbox"/>	Current (A)
<input type="checkbox"/>	Elapsed Time (s)
<input type="checkbox"/>	I Range (A)
<input type="checkbox"/>	Charge (C)
<input type="checkbox"/>	Sync ADC Input (V)
<input type="checkbox"/>	Forward I (A)
<input type="checkbox"/>	Reverse I (A)
<input type="checkbox"/>	Delta I (F-R) (A)
<input type="checkbox"/>	Delta I (R-F) (A)
<input type="checkbox"/>	Applied Potential (V)
<input type="checkbox"/>	Frequency (Hz)
<input type="checkbox"/>	SQRT of Frequency
<input type="checkbox"/>	1/SQRT of Frequency
<input type="checkbox"/>	Z  (ohms)
<input type="checkbox"/>	Zre (ohms)
<input type="checkbox"/>	Zim (ohms)
<input type="checkbox"/>	Phase of Z (deg)
<input type="checkbox"/>	Y  (S)
<input type="checkbox"/>	Yre (S)

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## 4 Versa Studio软件操作

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The screenshot shows a Microsoft Excel spreadsheet with the following data:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	Potential	Current (A)																
2	0.000614	1.55E-07																
3	0.00092	1.15E-06																
4	0.002147	2.01E-06																
5	0.003374	3.16E-06																
6	0.004294	4.15E-06																
7	0.005828	5.15E-06																
8	0.006748	6.14E-06																
9	0.007668	7.15E-06																
10	0.008588	8.15E-06																
11	0.009508	9.15E-06																
12	0.010428	1.02E-05																
13	0.011348	1.11E-05																
14	0.012268	1.21E-05																
15	0.013188	1.31E-05																
16	0.014108	1.41E-05																
17	0.015028	1.51E-05																
18	0.015948	1.61E-05																
19	0.016868	1.71E-05																
20	0.017788	1.81E-05																
21	0.018708	1.91E-05																
22	0.019628	2.01E-05																
23	0.020548	2.11E-05																
24	0.021468	2.21E-05																
25	0.022388	2.31E-05																
26	0.023308	2.41E-05																
27	0.024228	2.51E-05																

A callout box with a blue border and a white background points to the data table. The text inside the callout box reads: "粘贴到Excel或Origin等软件中" (Paste into Excel or Origin software).

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




"Potentiostatic"  
"Open Circuit"  
"Potentiodynamic"  
"Cyclic Polarization"  
"Linear Polarization"  
"Galvanostatic"  
"Galvanodynamic"  
"Cyclic Voltammetry (Single and Multiple Cycles)"  
"Staircase Cyclic Voltammetry (Single and Multiple Cycles)"  
"Linear Scan Voltammetry"  
"Staircase Linear Scan Voltammetry"  
"Chronoamperometry"  
"Chronopotentiometry"  
"Chronocoulometry"  
"Tafel"

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
#### 6.1.1 Open Circuit: 开路电位

Or

Equivalent Weight (g)	0
Mass (g)	0
Remeasure OC Per Action	Yes
Measured Open Circuit	0

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- Common
  - Linear Scan Voltammetry

Properties for Linear Scan Voltammetry

Step Properties	Value	Versus
Initial Potential (V)	-1	vs Ref
Final Potential (V)	1	vs Ref

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Scan Rate (V/s)	0.02
Total Points	1000

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

Annotations:

- 扫描电位范围 (Scan Potential Range)
- 扫描速率最大10V/s (Maximum Scan Rate 10V/s)
- 每次扫描固定取1000个数据点 (Fixed 1000 data points per scan)

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
#### 6.1.3 Cyclic Voltammetry (Single): 单圈循环伏安

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Experiment Properties

Actions to be Performed:

- Common
  - Cyclic Voltammetry

Properties for Cyclic Voltammetry

Endpoint Properties	Value	Versus	Vertex Hold	Acquire at Hold
Initial Potential (V)	-1	vs Ref		
Vertex Potential (V)	1	vs Ref	0	Yes
Final Potential (V)	-1	vs Ref		

Scan Properties	Value
Scan Rate (V/s)	1
Total Points	1999

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

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
#### 6.1.4 Cyclic Voltammetry (Multi Cycles): 多圈循环伏安

Experiment Properties

Actions to be Performed:

- Common
  - Cyclic Voltammetry (Multiple Cycles)

Properties for Cyclic Voltammetry (Multiple Cycles)

Endpoint Properties	Value	Versus	Vertex Hold	Acquire at Hold
Initial Potential (V)	0.5	vs Ref		
Vertex 1 Potential (V)	1	vs Ref	0	Yes
Vertex 2 Potential (V)	-1	vs Ref	0	Yes
Final Potential (V)	0	vs Ref		

Scan Properties	Value
Scan Rate (V/s)	0.1
Cycles	10
Points Per Cycle	2000
Total Points	22000

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto

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循环(Vertex 1→Vertex  
2→Vertex 1)次数

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
#### 6.1.5 Staircase Linear Scan Voltammetry: 阶跃线性扫描伏安

Experiment Properties

Actions to be Performed:

- Common
  - Staircase Linear Scan Voltammetry

Properties for Staircase Linear Scan Voltammetry

Endpoint Properties	Value	Versus	Limits	Direction	Value
Initial Potential (V)	0	vs Ref	None	≤	0
Final Potential (V)	1	vs Ref	None	≤	0

Scan Properties	Value
Step Height (mV)	10
Step Time (s)	1
Scan Rate (mV/s)	10
Total Points	101

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto

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
Experiment Properties

Actions to be Performed:

- Common
  - Staircase Cyclic Voltammetry

Properties for Staircase Cyclic Voltammetry				
Endpoint Properties	Value	Versus	Vertex Hold	Acquire at Hold
Initial Potential (V)	-1	vs Ref		
Vertex Potential (V)	1	vs Ref	5	Yes
Final Potential (V)	-1	vs Ref		

Scan Properties	Value
Step Height (mV)	10
Step Time (s)	1
Scan Rate (mV/s)	10
Total Points	407

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto

More Options. More Solutions.


<http://www.par-solartron.com.cn>

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#### 6.1.7 Staircase Cyclic Voltammetry (Multiple Cycles) : 多圈阶跃循环伏安

Experiment Properties

Actions to be Performed:

- Common
  - Staircase Cyclic Voltammetry (Multiple Cycles)

Endpoint Properties		Value	Versus	Vertex Hold	Acquire at Hold
Initial Potential (V)	0.1	vs Ref			
Vertex 1 Potential (V)	1	vs Ref	0	Yes	
Vertex 2 Potential (V)	-1	vs Ref	0	Yes	
Final Potential (V)	0	vs Ref			

Scan Properties		Value
Step Height (mV)	10	
Step Time (s)	1	
Cycles	10	

Cell Properties		Value
Leave Cell ON	No	
Cell to Use	External	

Instrument F  
Current Ran  
Acquisition I  
Electromete  
E Filter  
I Filter  
Bandwidth I  
LCI Bandwid  
iR Compens

More Options. More Solutions.


<http://www.par-solartron.com.cn>

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#### 6.1.8 Multi-Vertex Scan : 多拐点阶跃扫描

Experiment Properties

Actions to be Performed:

- Common
  - Multi-Vertex Scan

Properties for Multi-Vertex Scan

Endpoint Properties	Value	Versus	Vertex Hold	Acquire at Hold
Initial Potential (V)	0	vs OC		
Vertex 1 Potential ...	1	vs OC	0	Yes
Not Used	0	vs OC	0	Yes
Final Potential (V)	0.1	vs OC		

Scan Properties	Value	Limits	Direction	Value
Step Height (mV)	10	None	≤	0
Step Time (s)	1	None	≤	0
Scan Rate (mV/s)	10			
Total Points	191			

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**Actions to be Performed:**

- Common
  - Chronoamperometry

**Properties for Chronoamperometry**

Step Properties	Value	Versus
Potential (V)	1	vs Ref

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	1000
Total Points	1000

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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- Chronocoulometry

Properties for Chronocoulometry

Step Properties	Value	Versus	Pre-Elect (s)
Potential (V)	1	vs Ref	10

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	100
Total Points	100

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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- Insert
- Up
- Down
- Remove
- Help
- Advanced

**Actions to be Performed:**

- Common
- Loop #1
  - Recurrent Potential Pulses
  - Recurrent Potential Pulses
  - Recurrent Potential Pulses
  - Recurrent Potential Pulses

**Properties for Recurrent Potential Pulses**

Pulse Properties	Value	Versus
Potential (V)	1	vs Ref

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	10
Total Points	10

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	Yes
Cell to Use	External

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Experiment Properties

Actions to be Performed:

- Common
- Loop #1
  - Recurrent Galvanic Pulses
  - Recurrent Galvanic Pulses
  - Recurrent Galvanic Pulses
  - Recurrent Galvanic Pulses

Properties for Recurrent Galvanic Pulses	
<b>Pulse Properties</b>	<b>Value</b>
Current (A)	0.01
<b>Scan Properties</b>	<b>Value</b>
Time Per Point (s)	1
Duration (s)	10
Total Points	10

Limits	Direction	Value
None	≤	0
None	≤	0

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Cell Properties	Value
Leave Cell ON	Yes
Cell to Use	External

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Actions to be Performed:

- Common
  - Fast Potential Pulses

Properties for Fast Potential Pulses

Pulse Properties	Value	Versus	Width (s)
Number of pulses	5		
Potential (V) 1	1	vs Ref	0.100
Potential (V) 2	-1	vs Ref	0.100
Potential (V) 3	0.5	vs Ref	0.050
Potential (V) 4	-0.5	vs Ref	0.050
Potential (V) 5	0.75	vs Ref	0.100

Instrument Properties	Value
Current Range	2mA
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Scan Properties	Value
Time Per Point (s)	0.025
Cycles	2
Total Points	32
Total Duration	0.8

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

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Actions to be Performed:

- Common
  - Fast Galvanic Pulses

Properties for Fast Galvanic Pulses

Pulse Properties	Value	Width (s)
Number of pulses	5	
Current (A) 1	0.0005	0.100
Current (A) 2	0.00075	0.100
Current (A) 3	0.001	0.050
Current (A) 4	0.00125	0.050
Current (A) 5	0.001	0.200

Scan Properties	Value
Time Per Point (s)	0.025
Cycles	2
Total Points	40
Total Duration	1

Limits	Direction	Value
None	≤	0
None	≤	0

Instrument Properties	Value
Current Range	2mA
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Navigation: Insert, Up, Down, Remove, Help, Advanced

类似于Fast Potential Pulses

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
#### 6.1.16 Square Wave Voltammetry : 方波伏安法

Experiment Properties

Actions to be Performed:  
Common  
Square Wave Voltammetry

Properties for Square Wave Voltammetry					
Endpoint Properties	Value	Versus	Limits	Direction	Value
Initial Potential (V)	0	vs Ref	None	≤	0
Final Potential (V)	-1	vs Ref	None	≤	0

Scan Properties	Value
Pulse Height (mV)	25
Step Height (mV)	10
Frequency (Hz)	100
Scan Rate (mV/s)	1000
Total Points	202

Instrument Properties	Value
Current Range	2mA
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto

**Scan Rate=Step Height\*Frequency**

More Options. More Solutions.

<http://www.par-solartron.com.cn>

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
#### 6.1.17 Differential Pulse Voltammetry : 差分脉冲伏安法

Experiment Properties

Actions to be Performed:  
Common  
Differential Pulse Voltammetry

Endpoint Properties			Value	Versus	Limits	Direction	Value
Initial Potential (V)	0	vs Ref			None	≤	0
Final Potential (V)	-1	vs Ref			None	≤	0

Scan Properties		Value
Pulse Height (mV)		25
Pulse Width (s)		1
Step Height (mV)		10
Step Width (s)		2
Scan Rate (mV/s)		5
Total Points		202

Instrument Properties		Value
Current Range		2mA
Electrometer Mode		Auto
E Filter		Auto
I Filter		Auto
Bandwidth Limit		Auto
LCI Bandwidth Limit		Auto

Scan Rate=Step Height/Step Width

More Options. More Solutions.


<http://www.par-solartron.com.cn>

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#### 6.1.18 Normal Pulse Voltammetry: 常规脉冲伏安法

Experiment Properties

Actions to be Performed:

- Common
  - Normal Pulse Voltammetry

Properties for Normal Pulse Voltammetry

Endpoint Properties	Value	Versus	Limits	Direction	Value
Initial Potential (V)	0	vs Ref	None	≤	0
Final Potential (V)	-1	vs Ref	None	≤	0

Scan Properties	Value
Pulse Width (s)	1
Step Height (mV)	10
Step Width (s)	2
Scan Rate (mV/s)	5
Total Points	202

Instrument Properties	Value
Current Range	2mA
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto

**Scan Rate=Step Height/Step Width**

More Options. More Solutions.


<http://www.par-solartron.com.cn>

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#### 6.1.19 Reverse Normal Pulse Voltammetry : 反常规脉冲伏安法

Experiment Properties

Actions to be Performed:

- Common
  - Reverse Normal Pulse Voltammetry

Endpoint Properties			Value	Versus	Limits	Direction	Value
Initial Potential (V)	0	vs Ref			None	≤	0
Final Potential (V)	-1	vs Ref			None	≤	0

Scan Properties	Value	Instrument Properties	Value
Pulse Width (s)	1	Current Range	2mA
Step Height (mV)	10	Electrometer Mode	Auto
Step Width (s)	2	E Filter	Auto
Scan Rate (mV/s)	5	I Filter	Auto
Total Points	202	Bandwidth Limit	Auto

Scan Rate=Step Height/Step Width

More Options. More Solutions.


<http://www.par-solartron.com.cn>

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#### 6.2.1 Corrsion Open Circuit----开路电位

The screenshot shows the 'Experiment Properties' window. On the left is a navigation pane with 'Insert', 'Up', 'Down', 'Remove', 'Help', and 'Advanced' buttons. The main area is titled 'Actions to be Performed:' and lists 'Corrosion Open Circuit', 'Measure Open Circuit', and 'Linear Scan Voltammetry'. A callout box points to the 'Drift Rate' property in the 'Scan Properties' table.

Properties for Corrosion Open Circuit	
Scan Properties	Value
Time Per Point (s)	0.2
Duration (s)	60
Drift Rate (mV/min)	10
Total Points	300

Instrument Properties	Value
Current Range	2mA
Acquisition Mode	Auto
Electrometer Mode	Auto
E Resolution	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Limits	Direction	Value
None	≤	0

Cell Properties	Value
Cell to Use	External

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
Common

- Linear Polarization Resistance (LPR)

Properties for Linear Polarization Resistance (LPR)

Step Properties	Value	Versus
Initial Potential (V)	-0.02	vs OCP
Final Potential (V)	0.02	vs OCP

Limits	Direction	Value
None	≤	0
None	≤	0

Scan Properties	Value
Step Height (mV)	0.1
Step Time (s)	0.6
Scan Rate (mV/s)	0.1666
Total Points	401

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

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- Common
  - Tafel

**Properties for Tafel**

Endpoint Properties	Value	Versus
Initial Potential (V)	-0.25	vs OCP
Final Potential (V)	0.25	vs OCP

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Step Height (mV)	0.5
Step Time (s)	3
Scan Rate (mV/s)	0.1666
Total Points	1001

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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- Common
  - Potentiodynamic

**Properties for Potentiodynamic**

Endpoint Properties	Value	Versus
Initial Potential (V)	-0.25	vs OC
Final Potential (V)	1.6	vs Ref

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Step Height (mV)	1
Step Time (s)	6
Scan Rate (mV/s)	0.1666
Total Points	1851

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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Actions to be Performed:

- Common
  - Cyclic Polarization

Properties for Cyclic Polarization

Endpoint Properties	Value	Versus	Threshold Property	Value
Initial Potential (V)	-0.1	vs OC	Threshold	Enabled
Vertex Potential (V)	1.2	vs Ref	Start Level (V)	0.2
Final Potential (V)	0	vs OC	Threshold (mA)	1

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Step Height (mV)	1
Step Time (s)	1
Scan Rate (mV/s)	1
Total Points	2501

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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Common

Galvanic Corrosion

Properties for Galvanic Corrosion


Scan Properties	Value
Time Per Point (s)	0.5
Duration (s)	200
Total Points	400

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Cell to Use	External

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

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Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Step Height (mA)	0.1
Step Time (s)	10
Scan Rate (mA/s)	0.01
Total Points	11

Instrument Properties	Value
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Insert  
 Up  
 Down  
 Remove  
 Help  
 Advanced

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#### 6.2.10 Zero Resistance Ammeter (ZRA)—零阻计

Experiment Properties

Actions to be Performed:

Common

ZRA

Properties for ZRA


Scan Properties	Value
Time Per Point (s)	.1
Duration (s)	200
Total Points	2000

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Cell to Use	External

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Actions to be Performed: Common


Properties for EN

Scan Properties	Value
Time Per Point (s)	.1
Segment Duration (s)	200
Number of Segments	5
Total Points	10000

Limits	Direction	Value
None	≤	0
None	≤	0


Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Cell Properties	Value
Cell to Use	External

More Options. more solutions.

#### Galvanic Corrosion Cell Connection

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
- Common
  - Split LPR

**Properties for Split LPR**

Endpoint Properties	Value	Versus	Limits	Direction	Value	Cell Properties	Value
Cathodic Max (V)	-0.02	vs OC	None	≤	0	Leave Cell ON	No
Anodic Max (V)	0.02	vs OC	None	≤	0	Cell to Use	External

Scan Properties	Value
Step Height (mV)	0.1
Step Time (s)	0.6
Rest Duration (s)	200
Rest Drift Rate (mV/min)	0
Scan Rate (mV/s)	0.16666666
Total Points	401

Instrument Properties	Value
Current Range	Auto
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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Actions to be Performed:

- Common
  - Galvanic Control LPR

Properties for Galvanic Control LPR

Step Properties	Value
Initial Current (mA)	0
Final Current (mA)	1

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Step Height (mA)	0.1
Step Time (s)	10
Scan Rate (mA/s)	0.01
Total Points	11

Instrument Properties	Value
Acquisition Mode	Auto
Electrometer Mode	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

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
Insert  
Up  
Down  
Remove  
Help  
Advanced

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#### 6.3.1 Open Circuit----开路电位

**Experiment Properties**

**Actions to be Performed:**

- Common
  - Open Circuit
  - Measure Open Circuit
  - Linear Scan Voltammetry

**Properties for Open Circuit**

Scan Properties	Value
Time Per Point (s)	0.2
Duration (s)	60
Drift Rate (mV/min)	10
Total Points	300

Instrument Properties	Value
Current Range	2mA
Acquisition Mode	Auto
Electrometer Mode	Differential
E Resolution	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Limits	Direction	Value
None	≤	0

Cell Properties	Value
Cell to Use	External

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
Scan Properties		Value
Point Spacing	Logarithmic	
Number of Points	30	
Points Per Decade	10	
Data Quality	1	
Measurement Delay (s)	0	

DC Properties		
Step or Scan	Step	
Potential (V)	0.1	vs OC

Instrument Properties		Value
Current Range	Auto	
Electrometer Mode	Single Ended	
Bandwidth Limit	Auto	
LCI Bandwidth Limit	Auto	

Cell Properties		Value
Leave Cell ON	No	
Cell to Use	External	

Frequency List (Hz)	
1000000.000000	
794328.234724	
630957.344480	
501187.233627	
398107.170553	
316227.766017	
251188.643151	
199526.231497	
158489.319246	
125892.541179	
100000.000000	
79432.823472	
63095.734448	
50118.723362	

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Properties for Potentiostatic EIS

AC Properties	Value
Start Frequency (Hz)	100
End Frequency (Hz)	100
Amplitude (mV RMS)	10

DC Properties	Value	Versus
Step or Scan	Scan	
Initial Potential (V)	0	vs OC
Final Potential (V)	0.2	vs OC

Scan Properties	Value
Point Spacing	Linear
Number of Points	20
Points Per Decade	10
Data Quality	1
Measurement Delay (s)	300
Scan Rate (V/s)	0.01

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Single Ended
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

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Experiment Properties

Actions to be Performed:

- Common
  - Potentiostatic EIS
    - Loop #1
      - Potentiostatic EIS

Properties for Potentiostatic EIS

AC Properties	Value
Start Frequency (Hz)	10000
End Frequency (Hz)	1
Amplitude (mV RMS)	10

DC Properties	Value	Versus
Step or Scan	Step	
Potential (V)	0.02	vs Previc

Scan Properties	Value
Point Spacing	Logarithmic
Number of Points	30
Points Per Decade	10
Data Quality	1
Measurement Delay (s)	0

Frequency List (Hz)
10000.000000
7943.282347
6309.573445

Control Panel: Insert, Up, Down, Remove, Help, Advanced

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
Figure 5. Electrochemical impedance versus potential diagrams measured for a bronze sample immersed in artificial rainwater at 25°C.

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Experiment Properties

Actions to be Performed: Mott-Schottky

Properties for Mott-Schottky		
AC Properties	Value	
Start Frequency (Hz)	1000	
End Frequency (Hz)	1000	
Amplitude (mv RMS)	10	

DC Properties	Value	Versus
Initial Potential (V)	0	vs Ref
Final Potential (V)	1	vs Ref
Number of Steps	20	

Instrument Properties	Value
Current Range	Auto
Electrometer Mode	Single Ended
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Frequency List (Hz)	
1000.000000	

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
Experiment Properties

Actions to be Performed:

- Energy Open Circuit
- Measure Open Circuit
- Linear Scan Voltammetry

Properties for Energy Open Circuit

Scan Properties	Value
Time Per Point (s)	0.2
Duration (s)	60
Drift Rate (mV/min)	10
Delta Resolution	10
Delta E (mV)	0
Min. Tot	300
Max	300

Instrument Properties	Value
Current Range	2mA
Acquisition Mode	Auto
Electrometer Mode	Differential
E Resolution	Auto
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Limits	Direction	Value
None	≤	0

Cell Properties	Value
Cell to Use	External

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
#### 6.4.2 Constant Potential----恒电压充/放电

Experiment Properties - Data Acquired Monday, July 23, 2012 at 8:54:36 AM

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Properties for Constant Potential		
Step Properties	Value	Versus
Potential (V)	-0.05	vs OC
Accrue Q	Yes	

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	600
Delta Resolution	5
Delta I (mA)	0.5
Delta Q (mAh)	1
Min. Total Points	601
Max Total Points	3000

Properties	Value
Current(A)	2A
Potential(V)	
Capacity(mAh)	
Aux Input(V)	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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#### 6.4.3 Constant Current----恒电流充/放电

Experiment Properties - Data Acquired Monday, July 23, 2012 at 11:33:49 AM

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Instrument Properties		Value
Acquisition Mode		Average
Electrometer Mode		Auto 4/4
E Filter		Average
I Filter		None
Bandwidth Limit		Auto
LCI Bandwidth Limit		Auto

Scan Properties		Value
Time Per Point (s)		10
Duration (s)		1000
Delta Resolution		100
Delta E (mV)		1
Delta Q (mAh)		0
Min. Total Points		101
Max Total Points		10000

Cell Properties		Value
Leave Cell ON		No
Cell to Use		External

Limits	Direction	Value
Potential(V)	≤	1.2
None	≤	0

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

#### 6.4.6 Current CCDPL (Cyclic Charge-Discharge with Potential Limitation)

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Experiment Properties - Data Acquired Thursday, July 19, 2012 at 5:25:42 PM

**Actions to be Performed:**

Common

- Loop #1
  - Constant Current
  - Constant Potential
  - Energy Open Circuit
  - Constant Current
  - Constant Potential
  - Energy Open Circuit

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Experiment Properties - Data Acquired Thursday, July 19, 2012 at 5:25:42 PM

**Actions to be Performed:**

- Common
  - Loop #1
    - Constant Current
    - Constant Potential
    - Energy Open Circuit
    - Constant Current
    - Constant Potential
    - Energy Open Circuit

**Properties for Constant Current**

Step Properties	Value
Current (mA)	-500

Instrument Properties	Value
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Scan Properties	Value
Time Per Point (s)	10
Duration (s)	10000
Delta Resolution	10
Delta E (mV)	1
Delta Q (mAh)	0
Min. Total Points	1001
Max Total Points	10000

Limits	Direction	Value
Potential(V)	≤	2.8
None	≤	0

Experiment Properties - Data Acquired Thursday, July 19, 2012 at 5:25:42 PM

**Actions to be Performed:**

- Common
  - Loop #1
    - Constant Current
    - Constant Potential
    - Energy Open Circuit
    - Constant Current
    - Constant Potential
    - Energy Open Circuit

**Properties for Constant Potential**

Step Properties	Value	Versus
Potential (V)	0	vs Previ...
Accrue Q	No	

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	10

Instrument Properties	Value
Current Range	2A
Acquisition Mode	Auto

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Experiment Properties - Data Acquired Saturday, July 14, 2012 at 11:14:26 AM

Actions to be Performed:		Properties for Constant Power				
Common		Step Properties	Value	Limits	Direction	Value
Loop #1	Constant Power	Power (W)	0.02	Potential(V)	≤	0.5
	Energy Open Circuit	Current	Discharging	None	≤	0
	Constant Current	Limit III (mA)	50			

Experiment Properties - Data Acquired Saturday, July 14, 2012 at 11:14:26 AM

Actions to be Performed:		Properties for Constant Current		Instrument Properties		Value	
Common		Step Properties	Value	Acquisition Mode	Auto		
Loop #1	Constant Current	Current (A)	0.05	Electrometer Mode	Auto		
	Energy Open Circuit			E Filter	Auto		
	Constant Power			I Filter	Auto		
	Energy Open Circuit			Bandwidth Limit	Auto		
				LCI Bandwidth Limit	Auto		

Limits	Direction	Value
Potential(V)	≥	2
None	≤	0

More Options. More Solutions.

<http://www.par-solartron.com.cn>

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#### 6.4.8 Resistance CCD (Cyclic Charge-Discharge) ----恒电阻充放电循环

Experiment Properties - Data Acquired Sunday, July 15, 2012 at 5:04:56 PM

**Actions to be Performed:**

Common

- [-] Loop #1
  - ... Constant Resistance
  - ... Energy Open Circuit
  - ... Constant Current
  - ... Energy Open Circuit

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Experiment Properties - Data Acquired Sunday, July 15, 2012 at 5:04:56 PM

**Actions to be Performed:**

**Properties for Constant Resistance**

Step Properties	Value	Limits	Direction	Value
Resistance (ohms)	50	Potential(V)	≤	1
		None	≤	0

Experiment Properties - Data Acquired Sunday, July 15, 2012 at 5:04:56 PM

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Instrument Properties	Value
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	10
Delta Resolution	10
Delta E (mV)	0
Delta Q (mAh)	0
Min. Total Points	10
Max Total Points	10

Limits	Direction	Value
Potential(V)	>	3
None	≤	0

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#### 6.4.9 Charge-Discharge---充放电循环

**Properties for Constant Potential**

Step Properties	Value	Versus	Limits	Direction	Value
Potential (V)	0	vs Previ...	Current(A)	≤	0.0001
Accrue Q	Yes		Capacity(mAh)	≥	100

**Properties for Constant Current**

Step Properties	Value	Scan Properties	Value
Current (mA)	-100	Time Per Point (s)	1
		Duration (s)	10
		Delta Resolution	10
		Delta E (mV)	0
		Delta Q (mAh)	0
		Min. Total Points	10
		Max Total Points	10

Limits	Direction	Value
Potential(V)	≤	2
None	≤	0

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Instrument Properties	Value	Cell Properties	Value
Current Range	2A	Leave Cell ON	No
Acquisition Mode	Auto	Cell to Use	External
Electrometer Mode	Differential		
E Filter	Auto		
I Filter	Auto		
Bandwidth Limit	Auto		
LCI Bandwidth Limit	Auto		

CC Limits	Direction	Value	CV Limits	Direction	Value
Potential(V)	≥	3.5	Capacity(mAh)	≥	1150
None	≤	0	None	≤	0

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Experiment Properties

Actions to be Performed:

- Common
  - Loop #1
    - Constant Current
    - Energy Open Circuit

Properties for Constant Current

Step Properties	Value
Current (mA)	10

Instrument Properties	Value
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	10
Delta Resolution	10
Delta E (mV)	0
Delta Q (mAh)	0
Min. Total Points	10
Max Total Points	10

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Limits	Direction	Value
None	≤	0
None	≤	0

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Common

- Loop #1
  - Constant Potential

**Properties for Constant Potential**

Step Properties	Value	Versus
Potential (V)	0	vs OC
Accrue Q	No	

Limits	Direction	Value
None	≤	0
None	≤	0

Cell Properties	Value
Leave Cell ON	No
Cell to Use	External

Scan Properties	Value
Time Per Point (s)	1
Duration (s)	10
Delta Resolution	10
Delta I (mA)	0
Delta Q (mAh)	0
Min. Total Points	10
Max Total Points	10

Instrument Properties	Value
Current Range	2A
Acquisition Mode	Auto
Electrometer Mode	Differential
E Filter	Auto
I Filter	Auto
Bandwidth Limit	Auto
LCI Bandwidth Limit	Auto
iR Compensation	Disabled

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