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Instrument → maintenance → Furnace Origin Position Adjustment

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Instrument → maintenance → Burner Origin Position Adjustment

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Instrument → Furnace Nozzle Position

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Check following items before you ignite.

- ① Check the rest of fuel gas and support gas (C₂H₂ cylinder should be replaced with new one, if it's primary pressure is below 0.5MPa).
- ② Check gas supply pressure (Required gas pressure during combustion; Fuel 0.09MPa, Support 0
- ③ Check the burner slot is not clogged (Check the inside of the burner head also).
- ④ Insert the burner head firmly into the burner fixing port.
- ⑤ Check the nebulizer is fixed to the chamber by the nebulizer retaining plate.
- ⑥ Check gas pipes, the drain pipe and the safety valve (backward) are connected to the chamber
- ⑦ Check that drain tank is filled with water.
- ⑧ Check the end of the drain tube in the drain container is not in water.
- ⑨ Set fuel gas flow rate as follows in Atomizer/Gas Flow Rate Setup page when you ignite.

AIR-C₂H₂ Flame : Fuel Gas 2.0L/min, Support Gas (AIR) 15.0L/min
 N₂O-C₂H₂ Flame: Fuel Gas 7.0L/min, Support Gas (N₂O) 11.0L/min
 Make sure that red feather of N₂O-C₂H₂ flame is more than 2 to 3mm

*Set the instrument to flame analysis mode by selecting element or opening a file before you
 *If the flame should not be extinguished even by pressing the EXTINGUISH button, turn off the main unit POWER switch to shut off the gas supply and close the gas main valve.

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WizAard

! You can check the safty device now. Do you wish to check the fuel gas pressure monitor?

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WizAard

! Close the main valve of the Fuel Gas nearest to the instrument, then press Purge button to discharge the Gas inside of the instrument. Make sure to turn on the ventilation to purge the discharged fuel gas from the room. If the pipe between the main valve and the insrument is long, take time and press the purge button several times to purge all gas.

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WizAard

! You can check the safty device now. Do you wish to check the support gas pressure monitor?

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